

Merge VNA

14.3

[IHE Integration Statement](#)

© Copyright Merge Healthcare Solutions Inc. 2023.

Licensed materials – Property of Merge Healthcare Solutions Inc.

The content of this document is confidential information of Merge Healthcare Solutions Inc. and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare and the Merge Healthcare logo are trademarks of Merge Healthcare Inc.

Microsoft, Windows, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both.

All other names are trademarks or registered trademarks of their respective companies.

U.S. GOVERNMENT RESTRICTED RIGHTS:

This product is a “Commercial Item” offered with “Restricted Rights.” The Government’s rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation (“FAR”) 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare Solutions Inc.

For application support or to report issues with user documentation, contact Customer Support:

☎ 1-877-741-5369 (North America)
+31.20.514.5073 (Europe, the Middle East and Africa)
1800 316 746 (Australia)

✉ MergeSupport@merative.com

Part	Date	Revision	Description
ECM-43230	February 2023	1.0	Initial version for 14.3

The latest version of this document can be found at <https://mergecustomerforce.com/mergeusercommunity/login>.

IHE Integration Statement

The following table outlines the IHE Integration Statement.

IHE Integration Statement		Date
		2023
Vendor	Product Name	Version
Merge Healthcare Incorporated	Merge VNA	14.3
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager / Image Archive	None
Presentation of Grouped Procedures	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Evidence Documents	Image Manager / Image Archive	None
Key Image Note	Image Manager / Image Archive	None
Simple Image and Numeric Report	Report Repository	None
	Enterprise Report Repository	None
	External Report Repository Access	None
Access to Radiology Information	Image Manager / Image Archive	None
	Report Repository	None
	External Report Repository Access	None
Nuclear Medicine	Image Manager / Image Archive	None
Cardiac Catheterization Workflow	Image Manager / Image Archive	Intermittently Connected Modality, Cardiac Catheterization
	Performed Procedure Step Manager	None

Echo Workflow	Image Manager / Image Archive	Echocardiography, Intermittently Connected Modality
	Performed Procedure Step Manager	None
Evidence Documents - Cardiology	Image Manager / Image Archive	None
Consistent Time (CT)	Time Client	None
Mammography	Image Manager / Image Archive	None
Stress Testing Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Image Document Source	Set of DICOM Instances
	Document Repository	None
Cross-Enterprise Document Sharing (XDS.b)	Document Registry	Patient Identity Feed (HL7 V2), Patient Identity Feed (HL7 V3), On-Demand Documents, Reference ID, Document Metadata Update
	Document Repository	None
	Integrated Document Source/ Repository	Document Replace
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Patient Identifier Cross-reference (PIX)	Patient Identity Cross-reference Manager	None
	Patient Identifier Cross-reference Consumer	PIX Update Notification
Patient Demographics Query (PDQ)	Patient Demographics Supplier	None
	Patient Demographics Consumer	None
Imaging Object Change Management (IOCM)	Image Manager/Archive	None
	Change Requester	None
Document Metadata Subscription (DSUB)	Document Metadata Notification Broker	Patient-Independent Subscription, Folder Subscription
	Document Metadata Publisher	None
Multi-Patient Queries (MPQ)	Document Registry	None
Web-based Image Capture (WIC)	Image Manager	None
Internet address for vendor's IHE information	https://www.merative.com/support/imaging-product-support-details/merge-healthcare-statements	

Links to Standards Conformance Statements for the Implementation		
HL7	Refer to Merge Communities -- http://www.merge.com/Support/Support-Portal.aspx	
DICOM	https://www.merative.com/support/imaging-product-support-details/merge-healthcare-statements	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org