

Merge Cardio™

12.4

IHE INTEGRATION STATEMENT

© Copyright Merge Healthcare Solutions Inc. 2023.

Licensed materials - Property of Merge Healthcare Solutions Inc..

The content of this document is confidential information of Merge Healthcare Solutions Inc. and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare and the Merge Healthcare logo are trademarks of Merge Healthcare Inc.

Microsoft, Windows, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both.

All other names are trademarks or registered trademarks of their respective companies.

U.S. GOVERNMENT RESTRICTED RIGHTS:

This product is a "Commercial Item" offered with "Restricted Rights." The Government's rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation ("FAR") 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare Solutions Inc.

Manufacturer:

Merge Healthcare Incorporated

900 Walnut Ridge Drive

Hartland, WI 53029 USA

For application support or to report issues with user documentation, contact Customer Support:

☎ 1-877-741-5369 (North America)

✉ mergesupport@merative.com

Part	Date	Revision	Description
VRCS-36729	August 2023	1.0	Release 12.4.

The latest version of this document can be found at <https://merge.my.site.com/mergeusercommunity>.

IHE Integration Statement		Date
		2023
Vendor	Product Name	Version
Merge Healthcare Solutions Inc.	Merge Cardio	12.4
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager / Image Archive	
	Performed Procedure Step Manager	
	Image Display	
	Evidence Creator	Creator Performed Procedure Step
Patient Information Reconciliation	Image Manager / Image Archive	
	Performed Procedure Step Manager	
Consistent Presentation of Images	Image Manager / Image Archive	
	Image Display	
	Evidence Creator	
	Print Composer	
Presentation of Grouped Procedures	Image Manager / Image Archive	
	Performed Procedure Step Manager	
Evidence Documents	Image Manager / Image Archive	
	Image Display	
Key Image Note	Image Manager / Image Archive	
	Image Display	
	Evidence Creator	
Simple Image and Numeric Report	Report Repository	
	Report Reader	
	Enterprise Report Repository	
	External Report Repository Access	

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager / Image Archive	Multiple Sources
	Report Repository	Multiple Sources
	Image Display	Multiple Sources
	Report Reader	Multiple Sources
	External Report Repository Access	
Nuclear Medicine	Image Manager / Image Archive	
Cardiac Catheterization Workflow	Image Manager / Image Archive	Intermittently Connected Modality Cardiac Cath
	Performed Procedure Step Manager	
Echo Workflow	Image Manager / Image Archive	Echocardiography Intermittently Connected Modality
	Performed Procedure Step Manager	
Evidence Documents - Cardiology	Image Manager / Image Archive	
Retrieve ECG for Display	Display	
Retrieve Information for Display	Display	
Consistent Time	Time Client	
Patient Synchronized Applications	Patient Context Participant	
Portable Data for Imaging	Portable Media Creator	
	Portable Media Importer	
Mammography	Image Manager / Image Archive	
Stress Testing Workflow	Image Manager / Image Archive	
	Performed Procedure Step Manager	
Internet address for vendor's IHE information: https://www.merative.com/support/imaging-product-support-details/merge-healthcare-statements		
Links to Standards Conformance Statements for the Implementation		
DICOM	https://www.merative.com/support/imaging-product-support-details/merge-healthcare-statements	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org